

# **P25USB P30USB photodetector modules**

## **Sens-Tech Limited**

6A Langley Business Centre  
Station Road, Langley, Berkshire, SL3 8DS, UK

tel: +44 (0)1753 214714

fax: +44 (0)1753 214715

email: [info@sens-tech.com](mailto:info@sens-tech.com)

**[www.sens-tech.com](http://www.sens-tech.com)**

registered in England 668759

an ISO 9001 registered company

© Sens-Tech Limited, 2010

the contents may not be copied or disclosed to any  
unauthorised third party without written permission

# **user manual**

**issue 2 september 2010**



## P25USB P30USB photodetector modules

		page
<b>1</b>	<b>precautions</b>	
	1.1 light exposure	2
	1.2 supply voltage	2
	1.3 maintenance	2
	1.4 cleaning	2
	1.5 environment	2
	1.6 magnetic fields	2
<b>2</b>	<b>installation</b>	
	2.1 system requirement	3
	2.2 software installation	3
	2.3 interfaces	3
	2.4 light collection	3
<b>3</b>	<b>start-up software</b>	
	3.1 selecting module	4
	3.2 control window	5
	3.3 counter/timer window	7
<b>4</b>	<b>programming</b>	
	4.1 command description	9
<b>5</b>	<b>troubleshooting</b>	12
<b>6</b>	<b>module design</b>	13
<b>7</b>	<b>photon counting techniques</b>	
	7.1 operating principles	14
<b>8</b>	<b>appendix</b>	
	appendix A - interface characteristics	17
	appendix B - environmental conditions	18
	appendix C - outline drawings	19
	appendix D - the ActiveX <sup>®</sup> control	20

## 1 precautions

### 1.1 light exposure

The photodetector module is extremely sensitive to light. It is recommended that the module is stored in the dark. Exposure to bright light, such as fluorescent lighting or daylight, should be minimised. After exposure to bright light, the dark current or dark count rate may take up to 24 hours to return to the expected low level.

**note: the photodetector module must not be switched on (using the <HV on> command) in bright light; this may result in permanent damage to the product.**

### 1.2 supply voltage

Ensure the +5 V supply voltage does not exceed +5.25 V. Exceeding +6.0 V on the module or reversing the supply, may result in permanent damage (**see appendix A, module interface characteristics**).

### 1.3 maintenance

Photodetector modules contain no user-maintainable components. Because of the High Voltages present inside the package, dis-assembly must not be attempted by the user. Photodetector modules must be returned to **Sens-Tech** for repair.

### 1.4 cleaning

The window of the photodetector module may be cleaned, using lens tissue and alcohol, but only while it is disconnected from its power supply. The window should be treated as a normal optical component. The housing may also be cleaned with alcohol and lens tissue. The use of other cleaning agents is not recommended. When cleaning, observe the precautions stated under **light exposure**.

### 1.5 environment

The photodetector module must not be exposed to levels outside those specified in **appendix B, environmental conditions**.

### 1.6 magnetic fields

The module housing is made of mu-metal, providing good shielding from transverse magnetic fields. For applications involving strong axial magnetic fields, please ask for assistance from our technical sales staff.

## 2 installation

### 2.1 system requirements

- 1 PC operating under Microsoft Windows® XP or Windows®7 with a CD ROM drive
- 2 +5V power supply

### 2.2 software installation


Software has been provided with the module to allow immediate use of the product. However, if custom user programmes are required then **section 4, programming** defines the necessary interface protocol for users to follow.

Insert the **counter/timer software compact disc** into the CD-ROM drive of your computer. The set-up program will install the counter/timer software on your PC together with a virtual COM port driver for the USB interface.

**note: installation set-up starts automatically on most computers. If it does not run the set-up programme on the CD.**

### 2.3 interfaces (see appendix A, module interface characteristics)

- 2.3.1 connect the module to a light-tight housing or sample chamber.
- 2.3.2 connect the USB lead to the computer.
- 2.3.3 connect the power leads to the +5 V power supply. (red = +5 V, black = 0 V)
- 2.3.4 if used, make the appropriate connections to the **trigger Input** (grey) and **user output** (violet).



- 2.3.5 switch on the +5 V power supply.

**note: the photodetector module must be installed in a light-tight housing or sample chamber before the selecting the <HV on> command; failure to do so may damage it.**

### 2.3.6 Driver installation


After installation of the software, when the first module is connected to a USB port, the PC will detect new hardware and install the two drivers. Follow the on-screen instructions to carry out this process.

This process will repeat the first time a module is plugged into each of the PC's remaining USB ports.



### 3.2 control window

The control window opens automatically when a results window is opened. It may also be opened by selecting **Window** and then **Show control window**.



Note: If multiple modules are used, there will only be one control window and this will control all the modules in unison. If you wish to control multiple modules individually you should open an instance of the program for each package.

In the **Control window** set up the following:

#### 3.2.1 external trigger

Select the external trigger mode as either **Off**, **Fixed** or **Continuous**. These operate as follows;

- Off                    When the start button is clicked counting begins immediately. Counting will then continue until the set number of readings have been taken or the STOP button is clicked.
- Fixed                When the start button is clicked counting will not begin until the first falling edge on the external trigger input (grey wire). Counting will then continue until the set number of readings have been taken or the STOP button is clicked.
- Continuous        When the start button is clicked readings will only be taken whilst the trigger input is at logic high. Readings begin within 50us after the trigger input is taken high. When the trigger input is taken low the current reading period will complete before the gate is closed. Counting will continue until the set number of readings have been taken or the STOP button is clicked.

#### 3.2.2 readings

Select the total number of readings you wish to take (32000 max).

#### 3.2.3 period

Select the gating period of the reading (10 ms to 300 s).

By clicking on the box that states the units of the period setting, the user can toggle between ms (10 - 2500) and s (1 - 300). When ms (10 - 2500) is used, the software reads counts from the hardware dependent on the period set. When s (1 - 300) is used, the software reads counts from the hardware every 1 s and then accumulates the results until the desired period is achieved.

### 3.2.4 high voltage

Select the **HV voltage** as **Off** or **On** or **Follow start**, noting the **precautions** stated in section 1. The module will not count until the **HV voltage** has been activated. When **Follow start** is selected the HV will switch on automatically when the start button is clicked. (see also **Delay Start** below).

### 3.2.5 user signal


Select the **User Signal** as **Off** or **On** or **Follow start**. This will then control the external output line (violet wire) switching it between 0 volts (off) and +5 volts (on) to permit control of external apparatus, such as a shutter mechanism (**see appendix A, module interface characteristics**). When **Follow start** is selected the user output goes to +5V when the start button is clicked. This is unaffected by the setting of the delay start.

### 3.2.6 start / stop

Selecting **Start / Stop** will either **Start** or **Stop** counting. On selecting **Start**, the software will function in accordance with the selection made in the **Control Window**.

### 3.2.7 delay start

Setting a value here allows the start of the counting to be delayed by up to 10,000ms from clicking the start button. This can be used, for example, when using **Follow start** to allow the HV to settle before counting begins. The HV takes approximately 600ms to settle, however a warm up time of at least 2 seconds is recommended.



### 3.3 counter / timer software window

The following facilities are available from the Main Window, both as drop down menus and as items on the toolbar.

#### file:

- Connect** selects the module to be used, as stated in **section 3.1**.
- Open** allows the user to select a previously saved set of results for analysis.
- Save** allows the user to save a set of results.  
note: If there are multiple results windows open, the relevant window must be selected by clicking on it.
- Save as** allows the user to change the save file name.
- Exit** exits the program.

#### edit:

- copy** copies highlighted data to the clipboard to enable it to be pasted into other software applications, such as a spreadsheet or word processor.
- select all** selects and highlights all of the data in the table.

#### window:

- math** plots two existing sets of data as a ratio on a point by point basis.
- statistics** allows previously recorded data to be analysed to find min, max, mean and standard deviation.
- show control window** displays the counter control window

**help:**

**technical support**      this takes you to the **Sens-Tech** web site, where a list of technical support contacts are available.

**about**                      this states the software version being used.

## 4 programming

The tables on the following pages list the commands available to operate the module. This information is helpful when writing your own software.

In addition, an ActiveX<sup>®</sup> control driver is provided on the installation CD, together with example programmes for use with Excel<sup>®</sup>, Lab View<sup>®</sup>, Visual Basic<sup>®</sup> and Delphi<sup>®</sup>. Details of this module are shown in **appendix D, the ActiveX<sup>®</sup> control**.

### 4.1 command description

All commands sent from the PC must consist of strings of ASCII characters terminated with CR and LF, except the stop command which is just a CR. As shown by the following table, commands consist of a single upper case letter which will, in some cases, be followed by a value. The value is sent as a single ASCII character (except for command V which uses two). For example a decimal value of 21 would be transmitted as a single byte with the value 14 (hex). Values are limited to the range 0 to 255.

Commands will either be acknowledged with a pair of ASCII characters (as shown in the table) or will initiate a flow of data back to the PC. Where the command is not recognized there will be no reply. If the characters following a valid command are not as expected the reply BC (bad command) will be sent to the PC. Replies to the PC are not terminated with CR and LF.

The commands S, C, B, L and E cause the module to transmit readings back to the PC. The readings consist of blocks of four bytes which hold the binary value of each reading. The first byte transmitted is the most significant byte. If the most significant bit of the first byte is a 1 an overflow has occurred in the counter. The maximum possible uncorrected count for each reading is 67,108,863.

All commands will stop any reading in progress, or stop any EEPROM memory transmission to the PC, and then perform the specified action.

An **ST (Start Message)** is sent to the PC at power-up and whenever a **Reset** occurs.

command	action	returns to PC
	Power up or reset.	The <b>ST (start)</b> message.
<b>S</b>	Starts <b>R</b> readings where each reading is timed over P*10ms.	At the end of each reading, the reading is sent to the PC in binary form as 4 bytes with the most significant byte sent first. If the count exceeds 67,108,863 during a reading that reading is terminated immediately and an error reading is sent to the PC. The software then continues with the next reading.  An error reading is distinguished by having the most significant bit of the first byte as a 1, otherwise it is 0. Thus, a valid reading has a maximum of 31 bits.
<b>C</b>	Starts continuous readings where each reading is timed over P*10ms. The readings are terminated by sending the <b>Stop</b> command ( <b>CR</b> ) or other command	
<b>L</b>	Allows continuous readings. Readings are initiated, where each reading is timed over P*10ms, whenever the <b>External Trigger</b> pin is 1; and readings are inhibited whenever the <b>External Trigger</b> pin is 0.  Sending the <b>Stop</b> command ( <b>CR</b> ) or other command terminates this mode of operation.	Following the <b>VA (Valid)</b> message, the readings are sent to the PC as for the <b>S</b> and <b>C</b> commands above.
<b>E</b>	Allows <b>R</b> readings, where each reading is timed over P*10ms. Each reading is initiated at the falling edge of the external Trigger pin. This may be repeated any number of times. Sending the <b>Stop</b> command ( <b>CR</b> ) or other command terminates this mode of operation.	
<b>M</b>	Starts <b>R</b> readings where each reading is timed over P*10ms. The readings are sequentially stored in the non-volatile EEPROM along with the number of readings taken, overwriting the currently stored readings. Overflowed readings are also stored.  If <b>R</b> is greater than 126, only 126 readings are taken and stored.	No readings are sent to the PC, but when all of the readings have been taken the <b>MF</b> message is sent to the PC.
<b>B</b>	Sends the readings held in the non-volatile memory to the PC. These readings remain in the non-volatile memory until overwritten by the consequences of an <b>M</b> command.	The readings held in the non-volatile memory.
<b>CR</b> (carriage return)	Cancels the current set of readings.	The <b>SP (Stopped)</b> message.
<b>R</b> followed by a character.	Sets the <b>Number of Readings</b> to the byte value of the character following <b>R</b> . <b>A 0</b> is seen as 1. The default <b>Number of Readings</b> is 1.	The <b>VA (Valid)</b> message.
<b>P</b> followed by a character.	Sets the <b>Periods</b> setting <b>P</b> to the number following <b>P</b> . <b>A 0</b> is changed to 1. This sets the counter gate interval in multiples of 10ms, eg. P=100 gives a one second gate time. The default <b>Periods</b> setting is 10.	The <b>VA (Valid)</b> message.

command	action	returns to PC
<b>O</b> followed by a character.	Turns on the <b>Output Signal 1</b> (+5V out) when the character following the <b>O</b> has a byte value of 1, or turns off (0V out) the Output Signal 1 when it has a byte value of 0. Byte values greater than 1 cause a <b>BC</b> message to be sent to the PC The default value is off.	The <b>VA (Valid)</b> message.
<b>V</b> followed by two characters	Turns off the Voltage when both of the characters following the <b>V</b> have byte values of zero, or turns on the Voltage if either one of the characters has a byte value that is not zero (port 3.7 is 0) The default value is off.	The <b>VA (Valid)</b> message.
<b>D</b>	Turns on the high voltage.	The <b>VA (Valid)</b> message.

### example

To initialise and take 20 readings in fixed count mode using a 100ms gating period

Apply power to counter

Replies:           S T           After a power on reset, two ASCII characters are sent to PC to indicate startup initialisation is complete

PC sends:         R x CR LF       (Four ASCII characters), where x represents an ASCII character with the value in the range 1 to 255 setting the number of readings to be taken. In this example we want the character with the decimal value 20 and should transmit a single byte with the value 14 (hex).

replies:         V A           Two ASCII characters acknowledging a valid command

PC sends:         P y CR LF       (Four ASCII characters). This sets the gating period. y represents an ASCII character with a value in the range 1 to 255 and sets the gating period in multiples of 10ms. In this example we want 100ms and so y would be the ASCII character with the decimal value of 10 (10x10ms=100ms).

replies:         V A

PC sends:         S CR LF       (Three ASCII characters) This starts the measuring process. No V A reply is sent to PC.

replies:         b1 b2..b80       (80 bytes) The counter will transmit back a string of bytes containing the reading in binary . Twenty readings will produce 80 bytes in total. When the first reading is complete it will be transmitted as four binary bytes as detailed above.  
Whilst this transmission proceeds the next reading begins.

## 5 troubleshooting

### 1 no counts

check all connections have been made correctly

check the +5 V power supply is on and correctly set

check the HV has been selected to **on**

### 2 count overflow or very high count rates

check for light leaks

### 3 software will not install

check that the computer is running Windows® XP or Windows® 7.

### 4 no readings can be taken

check count mode not selected to **external trigger** or **trigger** not available.

### 5 communication timeout

check the module is powered up and connected to the USB port. If a module is momentarily powered down or disconnected, when in use, it will be necessary to restart the software before it can be re-selected.

## 6 module design

**P25USB** photodetector modules include a 25 mm diameter photomultiplier tube with a 22mm active photocathode diameter. The **P30USB** photodetector modules includes of a 30 mm photomultiplier tube with a 25 mm active photocathode diameter.

The photomultiplier is powered by a very efficient, +HV, Cockroft Walton power supply. Low power dissipation ensures negligible heat generation within the module.

The high voltage is factory set, at the optimum operating voltage, to ensure stability of operation and stable photon counting performance.

The signal from the photomultiplier is connected to a high performance amplifier discriminator combination, which is factory preset to accept signal but to reject electrical noise. This is described in **section 7, photon counting techniques**.

Pulses from the output of the amplifier/discriminator are passed to a counter which is gated on/off by the microcontroller according to the gate period setting. The count is latched through to a parallel to serial converter before continuing to the next counting period. This ensures the dead time of the counter is kept to an absolute minimum (11.1 $\mu$ s gap between gate periods).

The microcontroller has a maximum count rate capability of 100 MHz after linearity correction. All counts are automatically corrected for dead-time losses, caused by pulse pile up. This is particularly important at high count rates.

The module also incorporates a watchdog circuit such that in the event of internal software failure, the module will automatically reset itself to its default settings.

The module housing is mumetal<sup>®</sup>, providing good shielding from transverse magnetic fields.

## 7 photon counting techniques


### 7.1 operating principles:

These photodetector modules use the photomultiplier in the photon counting mode. The following paragraphs summarise the salient features of this technique and provide typical performance data for these modules.

#### 7.1.1 photon counting

The following section provides information on **photon counting**, the intended role for these modules.

**Photon counting** is the most sensitive method for measuring weak or short-lived light emissions. The photocathode has a wavelength-dependent probability of releasing a photoelectron for each incident photon. This probability is known as the **quantum efficiency (QE)** and is normally expressed as a percentage. Typical spectral response characteristics for the photocathodes used in photon counting detectors are shown in **figure 1**.



**figure 1**  
photomultiplier light sensitivity as a function of wavelength described as quantum efficiency (QE)

Photoelectrons emitted by the photocathode are accelerated and focussed onto the first dynode of the electron multiplier section of the photomultiplier. On impact, each photoelectron liberates a number of secondary electrons which are, in turn, accelerated and focussed onto the second dynode. The process is repeated at each subsequent dynode and the secondary electrons from the last dynode are collected at the anode.

Individual charge pulses from the photomultiplier vary in magnitude because of the statistical nature of the gain process but, so long as they exceed the preset threshold of the discriminator, each will be counted as one pulse. Only small pulses, resulting from single electrons originating from intermediate dynodes, will fall below the threshold and, thus, rejected.

The voltage applied to the photomultiplier in the detector is factory set to provide the optimum gain for photon counting. At this operating point on the signal counts vs HV plateau characteristic, the slope is less than 0.2 % per volt and, hence, we have stable performance.

### 7.1.2 signal-to-noise

The photoelectric effect is a quantum mechanical process subject to fluctuations described by Poisson statistics.

A steady light source generating  $m$  photoelectrons per second for a period of  $t$  seconds will produce an average of  $mt$  photoelectrons with a standard deviation of  $(mt)^{1/2}$ . Hence signal-to-noise ratio is given by:

$$mt/(mt)^{1/2} = (mt)^{1/2}$$

In **photon counting** there is negligible additional noise generated in the gain process or in the amplifier/discriminator hence this expression is also valid for the complete Detector.

It should be noted that the signal-to-noise ratio depends on the number of counts, hence for low count rates an increased measurement period will always be advantageous.

### 7.1.3. responsivity

The responsivity of the Detector is the output count rate, less dark counts, for a given light power incident on the photomultiplier window:

$$\text{responsivity} = \text{QE}(\lambda) \cdot f \cdot \left( \frac{\lambda}{h \cdot c} \right)$$

where  $\text{QE}(\lambda)$  is the quantum efficiency of the photocathode at wavelength  $\lambda$ ,


$h$  is the Planck Constant,  $6.626 \times 10^{-34}$  Js,

$c$  is the speed of light,  $3 \times 10^8$  ms<sup>-1</sup>,

and  $f$  is the collection efficiency of the first dynode, typically 0.95 for the standard range of photomultipliers used in photon counting.

Hence, when responding to light of 400 nm wavelength, a module with a bialkali photocathode, with typical QE of 0.25, has a responsivity of  $4.78 \times 10^{17}$  counts W<sup>-1</sup>, so 1 fW of 400 nm light incident on the photomultiplier window will produce a mean output count rate of 478 s<sup>-1</sup>, excluding dark counts.


The temperature coefficient of quantum efficiency, and hence of responsivity, is shown in **figure 2**.



**figure 2**  
temperature coefficient of various photocathodes used for photon counting


### 7.1.4 background

Background or dark count is the unwanted signal produced by the photomultiplier tube in the absence of light stimulation. Dark counts in photon counting detectors are a function of photocathode type and temperature as shown in **figure 3**.



**figure 3**  
The beneficial effect of cooling a photomultiplier

Typical dynamic range is shown in **figure 4**.



**figure 4**  
extension of the dynamic range by subtracting dark counts

Over the majority of their useful operating range the responsivity of Photon Counting Detectors is linear and will produce an output count rate which is proportional to input light power, as shown in **figure 4**.

At low count rates the effect of dark counts becomes significant causing the curves of figure 4 to flatten. Some improvement can be achieved by cooling, or at least preventing the module from being heated by other apparatus. Dark count subtraction can be used down to the point at which signal is less than about 1/10 of background. At lower signal rates the statistical uncertainty associated with the subtraction process becomes excessive.

In background subtraction, the signal count rate,  $S$ , is deduced as the difference between the measured rates of  $(S + B)$  for signal and background and  $(B)$  for background alone. It should therefore be noted that the signal-to-noise ratio for the deduced signal rate will be  $S/(S + 2B)^{1/2}$ .

At high count rates, the dynamic range is limited by electronic dead-time effects caused by pulse pile-up. Pulses that arrive whilst the discriminator is busy are not counted. Correction for dead-time is automatically applied in the software, in the form:

$$N = n/(1-nt)$$

where  $N$  is the corrected count rate,  $n$  is the measured count rate and  $t$  is the count rate correction factor.

## appendix A

### module interface characteristics

#### absolute maximum ratings

	test conditions	min	typ	max	units
supply voltage				6.0	V
user output load		-1.0		+7.0	V

#### dc characteristics (T = 25 °C, supply = +4.75 V to +5.25 V)

	test conditions	min	typ	max	units
<b>trigger input (TTL)</b>	<b>grey lead</b>				
input low volts (max)				0.9	V
input high volts (min)		1.9			V
<b>user output</b>	<b>violet lead</b>				
output low volts (max)	output sinking 20 mA			0.5	V
output high volts (min)	output sourcing 80 $\mu$ A	2.4			V
	output sourcing 12 $\mu$ A	4.5			V

## appendix B

### environmental conditions

**1 temperature (operating)**

+5 to +55 °C

**2 sinusoidal vibration (non-operating)**

10 to 500 Hz for 10 cycles in each axis, at a peak amplitude of 2 g, 1 octave/min.

**3 random vibration (non-operating)**

20 Hz at 0.02 g<sup>2</sup>/Hz  
to 50 Hz at 0.02 g<sup>2</sup>/Hz  
falling to 500 Hz at 0.001 g<sup>2</sup>/Hz

**4 humidity (non-condensing)**


93 % at 30 °C

**5 pressure (operating)**

can withstand pressure reductions from 68 kPa (0.68 bar) to 100 kPa (1 bar).

## appendix C

outline drawings (mm)



## appendix D

### the activeX® control

This is the documentation for the **ActiveX® control**. It describes the function and parameters of the methods, properties and events in the **control**.

Several examples are provided on the CD ROM to help in software development. Specifically, we provide examples for Visual Basic®, Labview®, Excel® and Delphi®. The examples are intended as a guide to programming with the **ActiveX® control** and not as complete solutions.

The following example describes the correct sequence for using the **ActiveX® control**.

#### 1 first, open the COM port

```
object.Open(1)
```

#### 2 then setup the device

```
object.Continuous = false          sets non-continuous mode ie fixed number readings.
```

```
object.Triggered = false          disable external trigger input.
```

```
object.OutputSignal = false       sets the user output (violet wire) to 0 volts.
```

```
object.OutputVoltage = true       turns on the photomultiplier HV supply.
```

```
object.Period = 100               sets gate period to one second.
```

```
object.ReadingCount = 100        sets number of readings to be taken as 100
```

#### 3 call the Start method and handle the events generated

```
object.Start                       starts counting
```

#### 4 finish by turning off the HV and closing the COM port

```
object.OutputVoltage = false       turns off the photomultiplier HV supply.
```

```
object.Close
```

#### remarks

it is not unusual for the com port to be open for the entire duration of the program.

**properties**

<b>BufferInUse (read only)</b>	Gets the number of readings currently in the buffer that are awaiting processing
<b>Buffersize</b>	Gets or sets the number of slots available to the software FIFO buffer
<b>Continuous</b>	Gets or sets the continuous mode of data acquisition. Setting this property to <b>TRUE</b> causes data acquisition to be continuous. Setting this value to false causes data acquisition to stop after the number of readings specified in the <b>Reading Count</b> property
<b>Output Signal</b>	Gets or sets the state of the output signal. Setting this property to <b>TRUE</b> causes the output signal set to +5 volts.
<b>Output Voltage</b>	Gets or sets the state of the photomultiplier high voltage supply. Setting this property to <b>TRUE</b> causes the high voltage to be switched on. (caution; see precautions in section 1.1)
<b>Period</b>	Gets or sets the counter gate period in multiples of 10ms.
<b>ReadingCount</b>	Gets or sets the number of readings to be taken in non-continuous acquisition mode
<b>Triggered</b>	Gets or sets the triggered mode of data acquisition. Setting this property to <b>TRUE</b> enables

**methods**

<b>Close</b>	Closes the <b>COM</b> port used by the ActiveX® control.
<b>Open</b>	Opens a <b>COM</b> port for use by the ActiveX® control.
<b>Start</b>	Begins counting.
<b>Stop</b>	Ends counting.

**events**

<b>Result</b>	Occurs when a new reading is available for processing
<b>OnBufferOverrun</b>	Occurs when the FIFO buffer attempts to grow to a size greater than the one specified in the <b>BufferSize</b> property.